

Pillion Rider Vario Footrest System



The Vario footrest system (joint) will be mounted to the motorcycle like the original footpeg. The Vario plate (adapter) will be mounted to the joint and the footpeg will be mounted to the Vario plate. You can only use our Super Grip footrest, our Racing footrest or our Multi Grip footrest. You cannot use the original footrest. You can fix the Vario plate every 45 degrees. For a better understanding you can see the complete footpeg lowering on pictures 2, 3 and 4 (all three parts must be ordered separately!). The joint is silver anodized.

Art-Nr.: 051 02 01 10

47,90 EUR

Pillion Rider Vario Footrest System (for left and right)

Additional pictures:

